



# SKIN HYDROPOWER®

The quintessence of  
dermal biostimulation  
and hydration



italfarmacia  
LABORATORI

Always first for your  
well-being!

# SKIN HYDROPOWER®

Dermal biostimulant with intense hydrating action

**HIGH CONCENTRATION OF MEDIUM  
AND HIGH MOLECULAR WEIGHT  
HYALURONIC ACIDE**

**“TARGETED” AMINO ACIDS**

# SKIN HYDROPOWER®

**HYALURONIC ACID** (medium and high molecular weight), able to significantly increase the level of skin hydration.

**GLYCINE**: a fundamental component of collagen; promotes elasticity, repair and firmness.

**PROLINE**: stabilizes collagen helices; contributes to mechanical strength and regeneration.

**VALINE**: stimulates protein synthesis, promotes cell regeneration, supports skin structure.

**LYSINE**: stimulates collagen production, helps reduce acne and herpes.

**ALANINE**: supports energy metabolism and protein synthesis; useful during metabolic stress.

**LEUCINE**: stimulates protein synthesis, promotes repair and muscle growth; useful for skin during the recovery phase.

**CYSTEINE**: antioxidant, keratin precursor; protects against free radicals, strengthens skin, hair and nails.

# SKIN HYDROPOWER®

**SKIN HYDROPOWER®**, thanks to the high concentration of hyaluronic acid with two molecular weights and the presence of “targeted” amino acids, ensures deep, long-lasting hydration.

**SKIN HYDROPOWER®** stimulates the production of collagen, elastin and hyaluronic acid, improving skin tone, firmness and elasticity.

**SKIN HYDROPOWER®** deeply hydrates the skin and counteracts the cytotoxic action of free radicals on fibroblasts.

**SKIN HYDROPOWER®** significantly improves skin structure.

# SKIN HYDROPOWER®

**SKIN HYDROPOWER®** gives the skin a brighter, more hydrated and plumper appearance, preventing and counteracting the effects of ageing. The treatment ensures natural, harmonious results, while respecting facial expressiveness.

**SKIN HYDROPOWER®** ensures strong, hydrated skin.





## COMPOSITION AND PACKAGING

**Class III medical device, MDR certified.**

1 blister containing 1 sterile, single-use, pre-filled 2 mL syringe, with Luer Lock connector and sterile needles 30G-13 mm: Sodium hyaluronate 25 mg/mL, L-isoleucine, L-leucine, L-lysine hydrochloride, L-proline, L-valine, glycine, L-serine, L-alanine, L-cysteine, phosphate buffer system, WFI (Water for Injections).

## CHEMICAL-PHYSICAL DATA

pH: 6.15–6.85

Osmolarity: 325–525 mOsm/kg

Loss on drying: 4.44–5.44%

Hyaluronic acid: 25 mg/mL

Molecular weight of hyaluronic acid:

1.2 and 2.4 MDa

Amino acids: 1.33–1.43%

## MICROBIOLOGICAL DATA

Sterile

Bacterial endotoxins < 0.05 EU/mL

## PHYSICAL DATA

Viscosity (cP):

Viscosity at 0.1/s: 135.0

Viscosity at 1/s: 67.0

Viscosity at 10/s: 12.5


# SKIN HYDROPOWER®

## INDICATIONS


**SKIN HYDROPOWER®**  
is indicated for the prevention and  
treatment of facial skin damage.

## INSTRUCTIONS FOR USE

1 session every 7 days, 4 times



1 session every 15 days, 2 times



Maintenance: 1 monthly session





## Always first for your well-being!

Via di Tor Sapienza, 7 Tel. +39 06 45441800  
00155 Roma Fax +39 06 45441801  
[www.italfarmacia.com](http://www.italfarmacia.com) [info@italfarmacia.com](mailto:info@italfarmacia.com)



 Made in Italy

